

**tynor**  
Ergonomically yours



Orthopaedic Appliances  
& Fracture Aids

## Body Belts & Braces



**Abdominal Support 9"**  
A-01 S,M,L,XL,XXL,  
XXXL,XXXXL



**TT /Abdominal Belt 8"**  
A-03 S,M,L,XL,XXL,XXXL



**LS Belt Lumbopore**  
A-04 S,M,L,XL,XXL



**Lumbo Sacral Belt**  
A-05 S,M,L,XL,  
XXL, XXXL



**Contoured L.S. Support**  
A-07 S,M,L,XL,XXL,  
XXXL,XXXXL



**Rib Belt**  
A-09 S,M,L,XL,XXL



**Chest Binder**  
A-11 S,M,L,XL,XXL,  
XXXL,XXXXL



**Ash Brace  
(Hyper Extension Brace)**  
A-12 S,L



**Taylor's Brace  
Short/Long**  
A-13 CH,UN,SPL



**Lumbo Support Uni**  
A-15 UN, SPL



**Hernia Belt**  
A-16 S,M,L,XL,XXL



**Pregnancy Back Support**  
A-20 S,M,L,XL,XXL



**Lumbo Lacepull Brace**  
A-29 UN, SPL



**Lacepull L.S. Belt**  
A-30 S,M,L,XL,XXL,XXXL

## Cervical Aids



**Collar Soft (Firm Density)**  
B-01 CH,S,M,L



**C Collar Soft w/ Support**  
B-02 S,M,L,XL



**C Collar Hard Adjustable**  
B-03 S,M,L



**Cervical Orthosis  
(Philadelphia) Ethafoam**  
B-05 S,M,L,XL



**Cervical Collar Soft**  
B-07 S,M,L,XL



**Cervical Pillow Regular**  
B-08 UN



**Cervical Orthosis  
(Philadelphia) Plastazote**  
B-10 CH,S,M,L,XL



**Contoured Cervical Pillow**  
B-19 UN



**C Collar Hard w/ Chin**  
B-20 S,M,L

## Fracture Aids



**P A Sling (Tropical)**  
C-01 CH,S,M,L,XL,XXL



**U Shoulder Immobiliser**  
C-02 CH,UN,SPL



**E Shoulder Immobiliser**  
C-03 S,M,L,XL,XXL



**Clavicle Brace w/ Buckle**  
C-04 CH,S,M,L,XL



**Clavicle Brace w/ Velcro**  
C-05 CH,S,M,L,XL,XXL



**P A Sling (Baggy)**  
C-06 CH,S,M,L,XL,XXL



**Pouch Arm Sling Urbane**  
C-21 S,M,L,XL



**Cast Shoe (Rocker Sole)**  
C-08 S,M,L,XL



**Cast Cover (Leg)**  
C-16 UN



**Cast Cover (Arm)**  
C-19 UN

## Ankle Supports



**Ankle Binder**  
D-01 S,M,L,XL



**Ankle Brace**  
D-02 CH,S,M,L,XL



**Anklet (Pair)**  
D-03 S,M,L,XL



**F. Drop Splint (R/L)**  
D-17 CH,S,M,L,XL



**Anklet Comfeel (Pair)**  
D-25 S,M,L,XL,XXL



**Ankle Splint**  
D-26 UN



**Air Ankle Splint**  
D-40 UN



**Walker Boot**  
D-32 CH,S,M,L,XL



**PF Night Splint Derotation**  
D-41 S,M,L



**F. Drop Splint w/ liner (R/L)**  
D-43 S,M,L,XL



**Walker Boot (Short)**  
D-45 S,M,L,XL

## Allied Products



**Weight Cuff**  
H-01 to 03 ½kg / 1kg / 2kg



**Exercising Ball**  
H-05 OR / NR



**Ortho Stockinete 10 Mt.**  
I-09 to 12 5,7,5,10,12,5cm



**Scrotal Support**  
I-59 S,M,L,XL



**Exercising Gel Ball**  
H-10 UN



**Cool Pack**  
H-11 UN



**Hot & Cold Pack**  
I-44 UN



**Coccyx Cushion Seat**  
H-13 UN



**Back Rest**  
I-46 UN



**Restrainer**  
I-65 UN



**Heating Pad Ortho**  
I-73 UN

## Wrist & Forearm Products



**Elastic Wrist Splint (R/L)**  
E-01 S, M, L, XL



**Wrist Brace w/ Double Lock**  
E-05 S,M,L



**Wrist Brace w/ Thumb**  
E-06 UN



**Wrist & Forearm Splint(R/L)**  
E-03 CH,S,M,L,XL



**Tennis Elbow Support**  
E-10 S,M,L,XL



**Forearm Splint**  
E-30 UN



**Wrist Splint (Ambidextrous)**  
E-43 S,M,L,XL



**Wrist Splint w/ Thumb**  
E-44 S,M,L,XL



**Hand Resting Splint(R/L)**  
E-29 CH,S,M,L



**Elbow Support**  
E-11 S,M,L,XL



**ROM Elbow Brace**  
E-46 UN

# Neoprene Products



**K Support Hinged (Neo)**  
J-01 S,M,L,XL,XXL,XXXL



**W Brace w/ Thumb (Neo)**  
J-03 UN



**Wrist Wrap (Neo)**  
J-04 UN



**Knee Wrap (Neo)**  
J-05 UN, SPL



**Abs Wrap (Neo)**  
J-06 UN,SPL



**OA Knee Support (Neo)(R/L)**  
J-08 S,M,L,XL,XXL



**K Support Sportif (Neo)**  
J-09 S,M,L,XL,XXL



**Ankle Support (Neo)**  
J-12 UN



**Shoulder Support (Neo)**  
J-14 UN,SPL



**Knee Wrap Hinged (Neo)**  
J-15 CH,S,M,L,XL,XXL



**Ankle Wrap (Neo)**  
J-16 UN

# Knee Supports



**Knee Cap (Pair)**  
D-04 S,M,L,XL,XXL,XXXL



**K Cap Open Patella**  
D-05 S,M,L,XL,XXL



**K.Cap (w/ Rigid Hinge)**  
D-06 S,M,L,XL,XXL



**Knee Cap w/ Patellar Ring**  
D-07 S,M,L,XL



**Elastic Knee Support**  
D-08 S,M,L,XL,XXL



**Functional Knee Support**  
D-09 S,M,L,XL,XXL



**R.O.M. Knee Brace**  
D-10 UN



**Knee Immobilizer 19"**  
D-11 S,M,L,XL,XXL



**Knee Immobilizer 14"**  
D-13 CH,S,M,L,XL,XXL



**Thigh Support**  
D-14 S,M,L,XL



**Pattelar Support**  
D-19 UN



**Knee Cap Comfeel (Pair)**  
D-23 S,M,L,XL



**Knee Immobilizer 22"**  
D-28 S,M,L,XL,XXL



**Knee Immobilizer Adj.**  
D-42 UN,SPL



**Knee Protector**  
D-46 S,M,L,XL,XXL

# Compression Stockings



**C Stocking Mid Thigh (Pair)**  
Classic, I-15 S,M,L,XL,XXL



**C Stocking Below Knee (Pair)**  
Classic, I-16 S,M,L,XL,XXL



**Tubular Elastic Bandage 1 Mt.**  
I-27 to 32 4,5,6,3,7,5,8,8,10 cm



**Anti Embolism Stocking**  
(DVT) Classic I-66, 69 S,M,L,XL



**Medical C Stocking (Pair)**  
Knee High I-67, 68 S,M,L,XL



**Medical C Stocking (Pair)**  
Thigh High I-70, 71 S,M,L,XL

## Traction Kits



**Skin Traction Set (PUF Liner)**  
G-01 CH,AD



**Leg Traction Brace**  
G-05 S,M,L



**Pelvic Traction Kit with Weight Bag**  
G-24 S,M,L,XL



**C.T. Kit (Sleeping) with Weight Bag**  
G-26 UN



**C.T. Kit (Sitting) with Weight Bag**  
G-25 UN



**Traction Weight Bag**  
G-20 UN

## Finger Splints



**Finger Cot**  
F-02 S,M,L



**Finger Ext. Splint**  
F-03 S,M,L



**Frog Splint**  
F-04 S,M,L



**Mallet Finger Splint**  
F-05 UN




**Thumb Spica Splint**  
F-06 UN


## Walking Aids



**Walking Stick**  
L-07 SIL, BLK



**Walking Stick (L type)**  
L-08 GOLD/SIL/BLK



**Walking Stick Quadripod**  
L-12 UN



**Quadra Stick**  
L-36 SIL



**Walking Stick Tripod**  
L-32 UN



**Auxiliary Crutch**  
L-21 S,M,L



**Elbow Crutch Adj.**  
L-13 CH, UN




**Pod Walking Stick (15,19,22 & 33mm)**  
L-09 to L-11 & L-33 UN




**Walker Invalid's**  
Recp. - All Pods  
L-20 CH, UN



**Walker Invalid's - Front Wheels-AL**  
L-28 UN



**Walker Invalid's HD - Front Wheels**  
L-29 UN



**Walker Invalid's HD - All Pods**  
L-34 UN



**Commode Chair**  
L-35 UN

## OAC Products



**Contoured L.S. Belt**  
L-01 S,M,L,XL,XXL



**Lumbo Sacral Belt**  
L-02 S,M,L,XL,XXL



**Abdominal Belt**  
L-03 S,M,L,XL,XXL





**P A Sling (Oxypore)**  
L-04 S,M,L,XL



**Heating Pad Ortho**  
L-23 UN (Size-Standard)










# Bio-metrics vs Product-metrics

## Body Size vs Product Size

BODY PART	MEASUREMENT	PRODUCT CODES	SIZE	INCHES	CM
Chest, Abdomen	 <p>Measure circumference around the waist immediately below navel</p>	A01, A05, A03, A07, L01, A04, L03, L02	S	28-32	70-80
			M	32-36	80-90
			L	36-40	90-100
			XL	40-44	100-110
			XXL	44-48	110-120
		A20	S	34-38	86-97
			M	38-42	97-107
			L	42-46	107-117
			XL	46-50	117-127
			S	32-36	80-90
	 <p>Measure circumference around the waist 3" below navel</p>	G06, G24	M	36-40	90-100
			L	40-44	100-110
			XL	44-48	110-120
			XXL	48-52	120-130
			S	28-32	70-80
	 <p>Measure circumference around the waist 2" below navel</p>	A16, I 59	M	32-36	80-90
			L	36-40	90-100
			XL	40-44	100-110
		I 59	Universal	28-42	70-106
			Spl. Size	42-52	106-132
			Universal	28-44	71-112
	 <p>Measure circumference around the waist below navel</p>	A15, J06, A13	Spl Size	44-52	112-132
			Universal	28-40	70-100
		A29	Spl Size	40-52	100-130
S & M			28-36	70-90	
A30		L & XL	36-44	90-110	
		XXL & XXXL	44-52	110-130	
		Short	upto 5'- 6"	upto 165	
 <p>Measure height of the patient</p>		A12	Long	Above 5'- 6"	Above 165
	Universal		4.6-5.6 Ft.	140-171	
	A13	Spl Size	5.6-6.6 Ft.	171-201	
		S	28-32	70-80	
 <p>Measure circumference around the chest</p>	A11, A09, C03	M	32-36	80-90	
		L	36-40	90-100	
		XL	40-44	100-110	
		XXL	44-48	110-120	
		XXXL	48-52	120-130	
	C 02	CH	Below 28	Below 71	
		Universal	28-42	70-107	
	J 14	Spl Size	42-52	107-130	
		Universal	26-36	66-91.4	
		Spl Size	36-48	91.4-121.9	
Neck	 <p>Measure distance between edge of sternum and lower jaw with neck in the extension position</p>	B01,B02,B03,B07,B10,B05,B20	CH	2-2.5	5.08-6.35
			S	3-3.5	7.5-8.75
			M	3.5-4	8.75-10
			L	4-4.5	10-11.25
			XL	4.5-5	11.25-12.5
		B08, B19	Universal	3-4.5	7.5-11.25
Forearm and Clavicles	 <p>Measure circumference around the chest below underarms</p>	C01,C06,L04,C21	CH	19.7-23.6	50-60
			S	27.6-31.5	70-80
			M	31.5-35.4	80-90
			L	35.4-39.4	90-100
			XL	39.4-43.3	100-110
		C04, C05	XXL	43.3-47.2	110-120
			CH	Below28	Below71
			S	28-34	71-86
			M	34-40	86-101
			L	40-46	101-116
XL	46-52	116-132			
XXL	52-58	132-147			

# Bio-metrics vs Product-metrics

## Body Size vs Product Size

Ankle	 <p>Measure circumference approx 2 inches above the inner ankle joint</p>	D01, D25, D02, D03, L05, D32	CH	Below 7.2	Below 18			
			S	7.2-8.4	18-21			
			M	8.4-9.6	21-24			
			L	9.6-10.8	24-27			
			XL	10.8-12	27-30			
			XXL	12-13.2	30-33			
J12, D26, J16			Universal	7.2-12	18-30			
Knee	 <p>Measure circumference around mid thigh - approx. 6 inches above knee</p>	D04, D06, D08, D09, J01, J09, J08, I15, D14, D13, D11, D28, D07, D23, J15, D07, D27, L06	CH	Below 14.8	Below 37.5			
			S	14.8-17.2	37-43			
			M	17.2-19.6	43-49			
			L	19.6-22	49-55			
			XL	22-24.4	55-61			
			XXL	24.4-26.8	61-67			
	I 69, I 70, I 71			S	14.6-16.9	37-43		
	J 05			M	16.5-18.9	42-48		
	D 10, D 42			L	18.5-20.9	47-53		
				XL	20.5-22.8	52-68		
				Universal	14.8-22	37-55		
				Spl Size	22-29.2	55-73		
				Universal	14.8-22	37-55		
	Knee	 <p>Measure circumference around calf - approx. 6 inches below knee</p>	I 16, G 05	S	8.5-11	21.5-27.9		
				M	11-13.5	27.9-34.2		
L				13.5-16	34.2-40.6			
XL				16-18.5	40.6-46.9			
I 66, I 67, I 68				S	9.8-11.8	25-30		
				M	11.4-13.4	29-34		
			L	13-15	33-38			
			XL	14.6-16.5	37-42			
D19			Universal	12-20	30-50			
Wrist		 <p>Measure circumference evenly around wrist joint</p>	E05, E 29, E01	S	4.9-6.1	12.5-15.5		
				M	6.1-7.3	15.5-18.5		
				L	7.3-8.5	18.5-21.5		
	XL			8.5-9.7	21.5-24.5			
	J03, J04, F06, E06, I65			Universal	5-9.5	12.5-23.8		
	 <p>Measure the width at palm</p>	E29	CH	2.35-2.75	6-7			
			S	2.6-3	6.6-7.5			
			M	3-3.3	7.6-8.5			
L			3.4-3.7	8.6-9.5				
Elbow	 <p>Measure circumference around the elbow in extended position</p>	E 10, E11, L24	S	7-10	18-25			
			M	10-12	25-30			
			L	12-14	30-35			
			XL	14-16	35-40			
Fore-Arm	 <p>Measure the circumference at a distance 6" from the wrist along the arm</p>	E03, E44	CH	5.9-7.2	15-18.3			
			S	7.2-8.5	18.3-21.6			
			M	8.5-9.8	21.6-24.9			
			L	9.8-11.1	24.9-28.2			
			XL	11.1-12.4	28.2-31.5			
			E46			Universal	7.2-11.1	18.3-28.2
E30			Universal	7.2-12	18.3-30.5			
Finger	 <p>Length of Fingers</p>	F02, F03, F04	S	Below 2	Below 5			
			M	2-3	5 - 7.5			
			L	3-4	7.5 - 10			
			F03, F05			Universal	Below 2 - 4	Below 5 - 10
Foot	 <p>Measure shoe size</p>	C08, C22, D17, D41, D43, K01, K02, K04, K08, K09, K10, K11, K12, K13	<b>Size</b>	<b>EUROPEAN</b>	<b>AMERICAN</b>			
			CH	25-30	8-12			
			CH3	25-27	8 - 9.5			
			CH2	28 - 30	10.3 - 12			
			CH1	31-33	12.8 - 2			
			S	34 - 36	2.8 - 4.4			
			M	37- 39	5.3 - 6.8			
			L	40- 42	7.5 - 9.0			
			XL	43-45	9.8 - 11.5			
			K14			Universal	34-45	2.8-11.75
			K15			Child	28-33	10.25-2
						Adult	34-39	2.75-6.75

# Silicone & Foot Care Products



**Insole Full Silicone (Pair)**  
K-01 S,M,L,XL



**Heel Cushion Silicone (Pair)**  
K-02 S,M,L,XL



**Heel Cup Silicone (Pair)**  
K-09 S,M,L



**Metatarsal Pad Silicone (Pair)**  
K-08 UN



**Toe Separator Silicone (Pair)**  
K-04 S,M,L



**Medial Arch Orthosis (Pair)**  
K-10 S,M,L,XL K-11 CH1, CH2, CH3



**Insole Omni Silicone (Pair)**  
K-13 S,M,L,XL



**Bunion Splint**  
K-14 UN



**Arch Support (Pair)**  
K-15 CH,AD



**Heel Protector (Pair)**  
K-17 UN



**Insole Gel Female (Pair)**  
K-16 S,M,L

# Child Range



**Collar Soft (Firm Density)**  
B-01 CH



**Cervical Orthosis Philadelphia**  
B-10 CH



**Pouch Arm Sling (Tropical)**  
C-01 CH



**Universal Shoulder Immobiliser**  
C-02 CH



**Clavical Brace with Buckle**  
C-04 CH



**Clavical Brace with Velcro**  
C-05 CH



**Pouch Arm Sling (Baggy)**  
C-06 CH



**Ankle Brace**  
D-02 CH



**Foot Drop Splint (R/L)**  
D-17 CH



**Walker Boot**  
D-32 CH



**Wrist & Forearm Splint**  
E-03 CH



**Hand Resting Splint**  
E-29 CH



**Skin Traction Set (PUF Liner)**  
G-01 CH



**Knee Wrap Hinged (Neo)**  
J-15 CH



**Knee Immobilizer**  
D-44 CH



**Medial Arch Orthosis**  
K-11 CH



**Arch Support**  
K-15 CH

ISO  
9001-2008



W.H.O.  
G.M.P.

ISO  
13485



US FDA  
registered